

## Product datasheet for **TA813874S**

### D-Dimer Mouse Monoclonal Antibody [Clone ID: OTI4A1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI4A1
Applications:	ELISA
Recommended Dilution:	ELISA 1:5000-10000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Native Human D-Dimer protein
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	180 kDa



[View online »](#)