

## **Product datasheet for TA814152S**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### TDP43 (TARDBP) Mouse Monoclonal Antibody [Clone ID: OTI4C1]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4C1
Applications: WB

Recommended Dilution: WB 1:500

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment of Human TARDBP (NP\_031401) produced in Ecoli.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Conjugation:** Unconjugated

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 44.7 kDa

**Gene Name:** TAR DNA binding protein

Database Link: NP 031401

Entrez Gene 230908 MouseEntrez Gene 23435 Human

Q13148

**Background:** The protein encoded by this gene is a transcriptional repressor that binds to chromosomally

integrated TAR DNA and represses HIV-1 transcription. In addition, this protein regulates alternate splicing of the CFTR gene. A similar pseudogene is present on chromosome 20.

[provided by RefSeq, Jul 2008]

Synonyms: ALS10; TDP-43





**Protein Families:** Transcription Factors

# **Product images:**





Figure A, Western blot analysis of overexpressed lysates(15ug per lane) from HEK293T cells transfected with empty plasmid ([PS100001], NC), human TARDBP plasmid ([RC210639], hTARDBP), mouse TARDBP plasmid ([MR206534], mTARDBP), rat TARDBP plasmid ([RR209466], rTARDBP) using anti-TARDBP antibody [TA814152](1:500). Figure B, Western blot analysis of the same samples as figure A with anti-DDK antibody ([TA180144], 1:1000).



Western blot analysis of extracts (50ug per lane) from 3 cell lines lysates by using anti-TARDBP monoclonal antibody([TA814152], 1:500).