

## Product datasheet for **TA814180S**

### **FABP3 Mouse Monoclonal Antibody [Clone ID: OTI11D7]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI11D7
Applications:	ELISA
Recommended Dilution:	ELISA 1:5000-10000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FABP3 (NP_004093) produced in EXPI-293F cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	14.9 kDa
Gene Name:	fatty acid binding protein 3
Database Link:	<a href="#">NP_004093</a> <a href="#">Entrez Gene 2170 Human</a> <a href="#">P05413</a>

**Background:** The intracellular fatty acid-binding proteins (FABPs) belongs to a multigene family. FABPs are divided into at least three distinct types, namely the hepatic-, intestinal- and cardiac-type. They form 14-15 kDa proteins and are thought to participate in the uptake, intracellular metabolism and/or transport of long-chain fatty acids. They may also be responsible in the modulation of cell growth and proliferation. Fatty acid-binding protein 3 gene contains four exons and its function is to arrest growth of mammary epithelial cells. This gene is a candidate tumor suppressor gene for human breast cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]



[View online »](#)

**Synonyms:** FABP11; H-FABP; M-FABP; MDGI; O-FABP

**Protein Pathways:** PPAR signaling pathway