

### **Product datasheet for TA814212**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CD21 (CR2) Mouse Monoclonal Antibody [Clone ID: OTI6C7]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI6C7

Applications: IHC

Recommended Dilution: IHC 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 208-587 of human

CD21 (NP 001868) produced in E.coli.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 113 kDa

**Gene Name:** complement component 3d receptor 2

Database Link: NP 001868

Entrez Gene 1380 Human

P20023

**Background:** This gene encodes a membrane protein, which functions as a receptor for Epstein-Barr virus

(EBV) binding on B and T lymphocytes. Genetic variations in this gene are associated with susceptibility to systemic lupus erythematosus type 9 (SLEB9). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep

2009]



#### CD21 (CR2) Mouse Monoclonal Antibody [Clone ID: OTI6C7] - TA814212


Synonyms: C3DR; CD21; CR; CVID7; SLEB9

**Protein Families:** Druggable Genome, Transmembrane


**Protein Pathways:** B cell receptor signaling pathway, Complement and coagulation cascades, Hematopoietic cell

lineage

# **Product images:**




Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-CD21 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, TA814212)



Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-CD21 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, TA814212)





Immunohistochemical staining of paraffinembedded Human spleen tissue within the normal limits using anti-CD21 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, TA814212)