

## Product datasheet for **TA814322**

### **IgG4 Mouse Monoclonal Antibody [Clone ID: OTI8D1]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI8D1
Applications:	ELISA
Recommended Dilution:	ELISA 1:1000-5000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Native Human IgG4 protein
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	150KD



[View online »](#)