

Product datasheet for TA814383

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

NGAL (LCN2) Mouse Monoclonal Antibody [Clone ID: OTI2G5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2G5
Applications: ELISA

Recommended Dilution: ELISA 1:5000-10000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 21-198 of human NGAL

(NP_005555) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Predicted Protein Size: 22.6 kDa

Gene Name: lipocalin 2

Database Link: NP 005555

Entrez Gene 3934 Human

P80188



NGAL (LCN2) Mouse Monoclonal Antibody [Clone ID: OTI2G5] - TA814383

Background: This gene encodes a protein that belongs to the lipocalin family. Members of this family

transport small hydrophobic molecules such as lipids, steroid hormones and retinoids. The protein encoded by this gene is a neutrophil gelatinase-associated lipocalin and plays a role in innate immunity by limiting bacterial growth as a result of sequestering iron-containing siderophores. The presence of this protein in blood and urine is an early biomarker of acute kidney injury. This protein is thought to be be involved in multiple cellular processes,

including maintenance of skin homeostasis, and suppression of invasiveness and metastasis. Mice lacking this gene are more susceptible to bacterial infection than wild type mice.

[provided by RefSeq, Sep 2015]

Synonyms: 24p3; MSFI; NGAL; p25

Protein Families: Secreted Protein