

## Product datasheet for **TA814427S**

### GFAP Mouse Monoclonal Antibody [Clone ID: OTI5H5]

#### Product data:

|                         |  |
|-------------------------|--|
| Product Type:           | Primary Antibodies   |
| Clone Name:             | OTI5H5   |
| Applications:           | ELISA  |
| Recommended Dilution:   | ELISA 1:5000-10000   |
| Reactivity:             | Human  |
| Host:                   | Mouse  |
| Isotype:                | IgG1   |
| Clonality:              | Monoclonal   |
| Immunogen:              | Human recombinant protein fragment corresponding to amino acids 292-432 of human GFAP (NP_002046) produced in E.coli.  |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.   |
| Concentration:          | 1 mg/ml  |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)  |
| Conjugation:            | Unconjugated   |
| Predicted Protein Size: | 49.9 kDa   |
| Gene Name:              | glial fibrillary acidic protein  |
| Database Link:          | <a href="#">NP_002046</a><br><a href="#">Entrez Gene 2670 Human</a><br><a href="#">P14136</a>  |
| Background:             | This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008] |
| Synonyms:               | ALXDRD   |
| Protein Families:       | ES Cell Differentiation/IPS  |



[View online »](#)