

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA814448

SARS-CoV-2 N Protein Mouse Monoclonal Antibody [Clone ID: OTI9F7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI9F7
Applications:	ELISA
Recommended Dilution:	ELISA 1:5000-10000
Reactivity:	SARS-CoV-2
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Full length recombinant SARS-Cov-2 N protein produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	47kDa
Database Link:	Entrez Gene 43740575 SARS-CoV-2
Background:	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive- sense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). Virus particles include the RNA genetic material and structural proteins needed for invasion of host cells. Once inside the cell the infecting RNA is used to encode structural proteins that make up virus particles, nonstructural proteins that direct virus assembly, transcription, replication and host control and accessory proteins whose function has not been determined.~ The structural proteins of SARS-CoV-2 include the envelope protein (E), spike or surface glycoprotein (S), membrane protein (M) and the nucleocapsid protein (N). The nucleocapsid phosphoprotein is a structural protein that binds to, protects the viral RNA genome and is involved in packaging the RNA into virus particles. The N protein has been suggested as an



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

antiviral drug target.