

## Product datasheet for **TA814477**

### **IgG4 Mouse Monoclonal Antibody [Clone ID: OTI1D10]**

#### **Product data:**

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Clone Name:             | OTI1D10   |
| Applications:           | ELISA   |
| Recommended Dilution:   | ELISA 1:5000-10000  |
| Reactivity:             | Human   |
| Host:                   | Mouse   |
| Isotype:                | IgG1  |
| Clonality:              | Monoclonal  |
| Immunogen:              | Native Human IgG4 protein   |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                                      |
| Concentration:          | 1 mg/ml   |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Conjugation:            | Unconjugated  |
| Predicted Protein Size: | 150KD   |



[View online »](#)