

## Product datasheet for **TA814510AM**

### **SARS-CoV-2 N Protein Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI6C3]**


#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI6C3
<b>Applications:</b>	ELISA, WB
<b>Recommended Dilution:</b>	WB 1:500-1000, ELISA 1:5000-10000
<b>Reactivity:</b>	SARS-CoV-2
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length recombinant SARS-CoV-2 N protein produced in HEK293.
<b>Formulation:</b>	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	0.5 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Conjugation:</b>	Biotin
<b>Storage:</b>	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Predicted Protein Size:</b>	47kDa
<b>Database Link:</b>	<a href="#">Entrez Gene 43740575 SARS-CoV-2</a>




[View online »](#)

Product images:



Western blot analysis of full length recombinant SARS-CoV-2 N protein (Cat# [TP790189], 0.02ug) by using anti-SARS-CoV-2 N protein antibody (Cat# [TA814510]).



Western blot analysis of full length recombinant SARS-Cov-2 N protein (Cat# [TP750180], 0.02ug) by using anti-SARS-CoV-2 N protein antibody (Cat# [TA814510]).