

Product datasheet for TA815031

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PAH Mouse Monoclonal Antibody [Clone ID: OTI4F6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4F6

Applications: WB

Recommended Dilution: WB 1:500

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PAH (NP_000268) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 51.8 kDa

Gene Name: phenylalanine hydroxylase

Database Link: NP 000268

Entrez Gene 18478 MouseEntrez Gene 5053 Human

P00439

Background: This gene encodes a member of the biopterin-dependent aromatic amino acid hydroxylase

protein family. The encoded phenylalanine hydroxylase enzyme hydroxylates phenylalanine to tyrosine and is the rate-limiting step in phenylalanine catabolism. Deficiency of this enzyme activity results in the autosomal recessive disorder phenylketonuria. [provided by

RefSeq, Aug 2017]





Synonyms: PH; PKU; PKU1

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Phenylalanine, tyrosine and tryptophan biosynthesis, Phenylalanine

metabolism

Product images:





Figure A, Western blot analysis of overexpressed lysates(15ug per lane) from HEK293T cells transfected with empty plasmid ([PS100001], lane 1), human PAH plasmid ([RC204694], lane 2), mouse PAH plasmid ([MR207240], lane 3), using anti-PAH antibody TA815031 (1:500). Figure B, Western blot analysis of the same samples as figure A with anti-DDK antibody ([TA180144], 1:1000)



Western blot analysis of extracts(50ug) from HepG2 cell lines lysates by using anti-PAH monoclonal antibody. (TA815031, 1:500)