

Product datasheet for TA815082S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Neuroligin 3 (NLGN3) Mouse Monoclonal Antibody [Clone ID: OTI6G10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI6G10

Applications: WB

Recommended Dilution: WB 1:500

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment of human NLGN3 (NP_061850) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 91.8 kDa

Gene Name: neuroligin 3

Database Link: NP 061850

Entrez Gene 245537 MouseEntrez Gene 54413 Human

O9NZ94

Background: This gene encodes a member of a family of neuronal cell surface proteins. Members of this

family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. Mutations in this gene may be associated with autism and Asperger syndrome. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Oct 2009]




Synonyms: HNL3

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs)

Product images:






Figure A, Western blot analysis of overexpressed lysates(15ug per lane) from HEK293T cells transfected with empty plasmid ([PS100001], lane 1), human NLGN1 plasmid ([RC207524], lane 2),human NLGN2 plasmid ([RC222544], lane 3), human NLGN3 plasmid ([RC207955], lane 4), human NLGN4X plasmid ([RC204671], lane 5), human NLGN4Y plasmid ([RC214651], lane 6), mouse NLGN3 plasmid ([MR222086], lane 7) using anti-NLGN3 antibody [TA815082] (1:500). Figure B, Western blot analysis of the same samples as figure A with anti-DDK antibody ([TA180144], 1:1000)

Western blot analysis of extracts (50ug) from mouse brain lysates by using anti-NLGN3 monoclonal antibody([TA815082], 1:500)