

Product datasheet for **TA890049S**

NUDT15 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human NUDT15
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3
Concentration:	1.64mg/ml
Purification:	Purified from the immunized serum by affinity chromatography (Protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18.4 kDa
Gene Name:	nudix hydrolase 15
Database Link:	NP_060753 Entrez Gene 55270 Human Q9NV35
Synonyms:	MTH2



[View online »](#)

Product images:



Western blot analysis of extracts (35ug) from different cell lines or tissues by using anti-NUDT15 rabbit polyclonal antibody ([TA890049]).