

Product datasheet for **TA890099S**


ANR52 (ANKRD52) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~2000
Reactivity:	Human, Rat, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human ANKRD52
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3
Concentration:	2.4 mg/ml
Purification:	Purified from the immunized serum by affinity chromatography (Protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	114.9 kDa
Gene Name:	ankyrin repeat domain 52
Database Link:	NP_775866 Entrez Gene 237615 MouseEntrez Gene 362811 RatEntrez Gene 283373 Human Q8NB46
Background:	Putative regulatory subunit of protein phosphatase 6 (PP6) that may be involved in the recognition of phosphoprotein substrates. [UniProtKB/Swiss-Prot Function]
Synonyms:	ANKRD33



[View online »](#)

Product images:


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ANKRD52 (Cat# [RC226317], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ANKRD52 antibody (Cat# [TA890099]). Positive lysates [LY430405] (100ug) and [LC430405] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from different cell lines and tissues by using anti-ANKRD52 rabbit polyclonal antibody.