

Product datasheet for **TA890112S**


DNMT3B Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human DNMT3B
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3
Concentration:	2.24 mg/ml
Purification:	Purified from the immunized serum by affinity chromatography (Protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	86 kDa
Gene Name:	DNA (cytosine-5-)-methyltransferase 3 beta
Database Link:	NP_787045 Entrez Gene 13436 MouseEntrez Gene 444985 RatEntrez Gene 1789 Human Q9UBC3
Synonyms:	ICF; ICF1; M.HsaIIIB
Protein Families:	Druggable Genome, Embryonic stem cells, Induced pluripotent stem cells, Stem cell - Pluripotency
Protein Pathways:	Cysteine and methionine metabolism, Metabolic pathways



[View online »](#)

Product images:


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY DNMT3B (Cat# [RC223206], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DNMT3B antibody (Cat# [TA890112]). Positive lysates [LY406229] (100ug) and [LC406229] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from different cell lines and tissues by using anti-DNMT3B rabbit polyclonal antibody.