

Product datasheet for **TA890119**


EPS15 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human EPS15
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3
Concentration:	4.57 mg/ml
Purification:	Purified from the immunized serum by affinity chromatography (Protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	98.5 kDa
Gene Name:	epidermal growth factor receptor pathway substrate 15
Database Link:	NP_001972 Entrez Gene 13858 MouseEntrez Gene 313474 RatEntrez Gene 2060 Human P42566
Synonyms:	AF-1P; AF1P; MLLT5
Protein Families:	Druggable Genome
Protein Pathways:	Endocytosis



[View online »](#)

Product images:


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY EPS15 (Cat# [RC212097], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EPS15 antibody (Cat# TA890119). Positive lysates [LY419614] (100ug) and [LC419614] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from different cell lines and tissues by using anti-EPS15 rabbit polyclonal antibody.