

## Product datasheet for **TA890121**

### **NMDAR2A (GRIN2A) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~2000
Reactivity:	Human, Rat, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human GRIN2A
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3
Concentration:	2.23 mg/ml
Purification:	Purified from the immunized serum by affinity chromatography (Protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	165.28 kDa
Gene Name:	glutamate ionotropic receptor NMDA type subunit 2A
Database Link:	<a href="#">NP_000824</a> <a href="#">Entrez Gene 14811 Mouse</a> <a href="#">Entrez Gene 24409 Rat</a> <a href="#">Entrez Gene 2903 Human</a> <a href="#">Q12879</a>

**Background:** This gene encodes a member of the glutamate-gated ion channel protein family. The encoded protein is an N-methyl-D-aspartate (NMDA) receptor subunit. NMDA receptors are both ligand-gated and voltage-dependent, and are involved in long-term potentiation, an activity-dependent increase in the efficiency of synaptic transmission thought to underlie certain kinds of memory and learning. These receptors are permeable to calcium ions, and activation results in a calcium influx into post-synaptic cells, which results in the activation of several signaling cascades. Disruption of this gene is associated with focal epilepsy and speech disorder with or without mental retardation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014]




[View online »](#)

**Synonyms:** EPND; FESD; GluN2A; LKS; NMDAR2A; NR2A


**Protein Families:** Druggable Genome, Ion Channels: Glutamate Receptors, Ion Channels: Sodium, Transmembrane

**Protein Pathways:** Alzheimer's disease, Amyotrophic lateral sclerosis (ALS), Calcium signaling pathway, Long-term potentiation, Neuroactive ligand-receptor interaction, Systemic lupus erythematosus

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY GRIN2A (Cat# [RC223136], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GRIN2A rabbit polyclonal antibody (Cat# TA890121).



Western blot analysis of extracts (35ug) from different cell lines and tissues by using anti-GRIN2A rabbit polyclonal antibody.