

# **Product datasheet for TA890226**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **SARS-CoV-2 N Protein Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type: Primary Antibodies

Applications: ELISA, ICC, IHC, WB

Recommended Dilution: ELISA 1:5000-1:10000, WB 1:2000, IHC 1:500-1:1000, ICC 1:500-1:1000

**Reactivity:** SARS-CoV-2

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**Immunogen:** Full length recombinant SARS-CoV-2 N protein produced in E.coli.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from the immunized serum by affinity chromatography (Protein A/G)

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.


**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 47kDa


**Database Link:** Entrez Gene 43740575 SARS-CoV-2



# **Product images:**




Western blot analysis of full length recombinant SARS-CoV-2 N protein(Cat# [TP790189], 0.02ug) by using anti-SARS-CoV-2 N protein antibody (Cat# TA890226). (1:2000)







Left: HEK293T transfected with plasmids overexpressing SARS-CoV-2 N gene (Cat# [VC102555]). Cells were fixed in 2% neutral buffer formalin and detected with N protein rabbit polyclonal antibody (Cat# TA890226). Polink1 anti rabbit HRP polymer (Cat# [D13-18]) was done for 30 min in room temperature followed by DAB. 40x image. Right: Untransfected HEK293T.





Screening for N-protein in paraffin-embedded on left mock transfected HEK 293T cells negative and on the right HEK 293T transfected with Myc-DDK-tagged ORF clone (Cat# [VC102555]) for SARS-CoV-2 nucleocapsid phosphoprotein showing positive stain. Immuno-cytochemical staining was done with rabbit polyclonal anti-SARS-CoV-2 N Protein (CAT# TA890226) at 1:200. No heat-induced epitope retrieval was used.





Screening for N-protein in paraffin-embedded on left mock transfected HEK 293T cells negative and on the right HEK 293T transfected with Myc-DDK-tagged ORF clone ([VC102555]) for SARS-CoV-2 nucleocapsid phosphoprotein showing positive stain. Immuno-cytochemical staining was done with rabbit polyclonal antiSARS-CoV-2 N Protein (TA890226) at 1:200 no heat-induced epitope retrieval was used.