

Product datasheet for **TA890227M**

SARS-CoV-2 Spike Protein Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB 1:2000, ELISA 1:5000-1:10000
Reactivity:	SARS-CoV-2
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein of SARS-CoV-2 spike RBD protein and virion surface domain of membrane protein (M-protein) produced in HEK293
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from the immunized serum by affinity chromatography (Protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34 kDa
Gene Name:	S Protein
Database Link:	Entrez Gene 43740568 SARS-CoV-2




[View online »](#)

Product images:



Western blot analysis of recombinant fusion protein of SARS-CoV-2 Spike RBD Protein and Virion Surface Domain of Membrane Protein ([TP701119], 0.02ug) by using anti-SARS-Cov-2 spike protein antibody ([TA890227]). (1:2000)



Western blot analysis of recombinant SARS-CoV-2 Spike RBD Protein ([TP701122], 0.1ug) by using anti-SARS-Cov-2 spike protein antibody ([TA890227]). (1:2000)