

Product datasheet for TA899002S

PMS2 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA, WB

Recommended Dilution: ELISA 1:5000-10000, WB 1:500

Reactivity: Human Rabbit Host: Isotype: lgG

Clonality: Polyclonal

Immunogen: This antibody is raised against a recombinant protein of PMS2 expressed in HEK293 cells

Formulation: Tissue culture buffer containing 1% BSA and 0.02% Sodium azide

Concentration: Lot dependent; please refer to CoA along with shipment

Purification: Purified from the immunized serum by affinity chromatography (Protein A/G)

Conjugation: Unconjugated

Store at -20°C as received. Storage:

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 95.8 kDa

Gene Name: PMS1 homolog 2, mismatch repair system component

Database Link: NP 000526

Entrez Gene 5395 Human

P54278

Background: This gene is one of the PMS2 gene family members found in clusters on chromosome 7. The

> product of this gene is involved in DNA mismatch repair. It forms a heterodimer with MLH1 and this complex interacts with other complexes bound to mismatched bases. Mutations in this gene are associated with hereditary nonpolyposis colorectal cancer, Turcot syndrome, and are a cause of supratentorial primitive neuroectodermal tumors. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2008]

Synonyms: HNPCC4; MLH4; MMRCS4; PMS2CL; PMSL2

Protein Families: Druggable Genome



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200


CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Protein Pathways: Mismatch repair

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PMS2 cDNA (Cat# [RC210880], Right lane) for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PMS2(Cat# [TA899002]). (1:500)