

Product datasheet for **TB400160**

C9orf95 (NMRK1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NMRK1 (NM_017881), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200160
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_017881</u> , <u>NP_060351</u>
Synonyms:	bA235O14.2; C9orf95; NRK1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_060351</u>
Locus ID:	54981
Cytogenetics:	9q21.13
Protein Pathways:	Nicotinate and nicotinamide metabolism



[View online »](#)