

Product datasheet for **TB400171**

DUSP23 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DUSP23 (NM_017823) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200171
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	DUSP23 Mouse Monoclonal Antibody [Clone ID: OTI3C10] (TA503944)
ACCN:	NM_017823 , NP_060293
Synonyms:	DUSP25; LDP-3; LDP3; MOSP; VHZ
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060293
Locus ID:	54935
Cytogenetics:	1q23.2
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)