

## Product datasheet for **TB400242**

### HRASLS3 (PLA2G16) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PLA2G16 (NM_007069), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200242
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HRASLS3 (PLA2G16) Mouse Monoclonal Antibody [Clone ID: OT11A5] (TA506908)
ACCN:	<a href="#">NM_007069</a> , <a href="#">NP_009000</a>
Synonyms:	AdPLA; H-REV107; H-REV107-1; HRASLS3; HREV107; HREV107-1; HREV107-3; HRSL3; PLA2G16; PLAAT-3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_009000</a>
Locus ID:	11145
Cytogenetics:	11q12.3-q13.1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)