

Product datasheet for **TB400439**

PPM1G CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PPM1G (NM_002707), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200439
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PPM1G Mouse Monoclonal Antibody [Clone ID: OT11C12] (TA502232)
ACCN:	NM_002707 , NP_002698
Synonyms:	MGC1675; MGC2870; PP2CG; PP2CGAMMA; PPP2CG
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002698
Locus ID:	5496
Cytogenetics:	2p23.3
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)