

## Product datasheet for **TB400440**

### **RPB2 (POLR2B) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of POLR2B (NM_000938) in HEK293T cells paraffin embedded controls for ICC/IHC staining
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC200440
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
<b>ACCN:</b>	<u>NM_000938</u> , <u>NP_000929</u>
<b>Synonyms:</b>	hRPB140; POL2RB; RPB2
<b>Storage:</b>	Room Temperature, or 2-8°C for long term storage
<b>Stability:</b>	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u>NP_000929</u>
<b>Locus ID:</b>	5431
<b>Cytogenetics:</b>	4q12
<b>Protein Families:</b>	Transcription Factors
<b>Protein Pathways:</b>	Huntington's disease, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase



[View online »](#)