

Product datasheet for **TB400680**

NME2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NME2 (NM_001018139), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200680
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NME2 Mouse Monoclonal Antibody [Clone ID: OTI2F9] (TA505574)
ACCN:	NM_001018139 , NP_001018149
Synonyms:	NDKB; NDPK-B; NDPKB; NM23-H2; NM23B; PUF
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001018149
Locus ID:	4831
Cytogenetics:	17q21.33
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Metabolic pathways, Purine metabolism, Pyrimidine metabolism



[View online »](#)