

Product datasheet for **TB400777**

hnRNP G (RBMX) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RBMX (NM_002139), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200777
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_002139 , NP_002130
Synonyms:	hnRNP-G; HNRNPG; HNRPG; MRXS11; RBMXP1; RBMXRT; RNMX
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002130
Locus ID:	27316
Cytogenetics:	Xq26.3
Protein Pathways:	Spliceosome



[View online »](#)