

Product datasheet for **TB401109**

HAGH CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HAGH (NM_005326), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201109
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_005326 , NP_005317
Synonyms:	GLO2; GLX2; GLXII; HAGH1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005317
Locus ID:	3029
Cytogenetics:	16p13.3
Protein Families:	Druggable Genome
Protein Pathways:	Pyruvate metabolism



[View online »](#)