

## Product datasheet for **TB401225**

### MMP7 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MMP7 (NM_002423) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201225
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	MMP7 Mouse Monoclonal Antibody [Clone ID: OTI3G7] (TA804842)
ACCN:	<a href="#">NM_002423</a> , <a href="#">NP_002414</a>
Synonyms:	MMP-7; MPLS1; PUMP-1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_002414</a>
Locus ID:	4316
Cytogenetics:	11q22.2
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Wnt signaling pathway



[View online »](#)