

Product datasheet for **TB401279**

beta Arrestin 1 (ARRB1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ARRB1 (NM_004041), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201279
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	beta Arrestin 1 (ARRB1) Mouse Monoclonal Antibody [Clone ID: OTI2B3] (TA812546)
ACCN:	NM_004041 , NP_004032
Synonyms:	ARB1; ARR1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004032
Locus ID:	408
Cytogenetics:	11q13.4
Protein Families:	Druggable Genome
Protein Pathways:	Chemokine signaling pathway, Endocytosis, MAPK signaling pathway



[View online »](#)