

## Product datasheet for **TB401294**

### TFIIF (GTF2F1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GTF2F1 (NM_002096) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201294
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	TFIIF (GTF2F1) Mouse Monoclonal Antibody [Clone ID: OTI4B10] (TA503604)
ACCN:	<u><a href="#">NM_002096</a></u> , <u><a href="#">NP_002087</a></u>
Synonyms:	BTF4; RAP74; TF2F1; TFIIF
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_002087</a></u>
Locus ID:	2962
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Basal transcription factors



[View online »](#)