

## Product datasheet for **TB401539**

### NDUFA5 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NDUFA5 (NM_005000) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201539
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	NDUFA5 Mouse Monoclonal Antibody [Clone ID: OTI1E8] (TA503881)
ACCN:	<u><a href="#">NM_005000</a></u> , <u><a href="#">NP_004991</a></u>
Synonyms:	B13; CI-13kB; CI-13KD-B; NUFM; UQOR13
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_004991</a></u>
Locus ID:	4698
Cytogenetics:	7q31.32
Protein Pathways:	Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease



[View online »](#)