

Product datasheet for **TB401544**

Integrin Linked Kinase (ILK) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ILK (NM_004517), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201544
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Integrin Linked Kinase (ILK) Mouse Monoclonal Antibody [Clone ID: OT12H3] (TA804883)
ACCN:	NM_004517 , NP_004508
Synonyms:	HEL-S-28; ILK-1; ILK-2; P59; p59ILK
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004508
Locus ID:	3611
Cytogenetics:	11p15.4
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Endometrial cancer, Focal adhesion, PPAR signaling pathway



[View online »](#)