

Product datasheet for **TB401569**

Annexin A1 (ANXA1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANXA1 (NM_000700) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201569
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Annexin A1 (ANXA1) Mouse Monoclonal Antibody [Clone ID: OTI3A8] (TA500967)
ACCN:	<u>NM_000700</u> , <u>NP_000691</u>
Synonyms:	ANX1; LPC1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_000691</u>
Locus ID:	301
Cytogenetics:	9q21.13
Protein Families:	Druggable Genome



[View online »](#)