

Product datasheet for **TB401661**

p27 KIP 1 (CDKN1B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CDKN1B (NM_004064) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201661
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_004064</u> , <u>NP_004055</u>
Synonyms:	CDKN4; KIP1; MEN1B; MEN4; P27KIP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_004055</u>
Locus ID:	1027
Cytogenetics:	12p13.1
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Chronic myeloid leukemia, ErbB signaling pathway, Pathways in cancer, Prostate cancer, Small cell lung cancer



[View online »](#)