

Product datasheet for **TB401681**

ARMER (ARL6IP1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ARL6IP1 (NM_015161) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201681
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_015161</u> , <u>NP_055976</u>
Synonyms:	AIP1; ARL6IP1; ARMER; SPG61
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_055976</u>
Locus ID:	23204
Cytogenetics:	16p12.3
Protein Families:	Transmembrane



[View online »](#)