

Product datasheet for **TB401713**

SNRP70 (SNRNP70) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SNRNP70 (NM_003089) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201713
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_003089 , NP_003080
Synonyms:	RNPU1Z; RPU1; Snp1; SNRP70; U1-70K; U1AP; U1RNP; U170K
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003080
Locus ID:	6625
Cytogenetics:	19q13.33
Protein Families:	Druggable Genome
Protein Pathways:	Spliceosome



[View online »](#)