

Product datasheet for **TB401836**

RPS4X CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RPS4X (NM_001007) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201836
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_001007 , NP_000998
Synonyms:	CCG2; DXS306; RPS4; S4; SCAR; SCR10
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000998
Locus ID:	6191
Cytogenetics:	Xq13.1
Protein Pathways:	Ribosome



[View online »](#)