

## Product datasheet for **TB401846**

### ARIH2 CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of ARIH2 (NM_006321) in HEK293T cells paraffin embedded controls for ICC/IHC staining  |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC201846  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)   |
| Target Detection Antibodies:          | ARIH2 Mouse Monoclonal Antibody [Clone ID: OTI3F9] (TA804094)   |
| ACCN:                                 | <u><a href="#">NM_006321</a></u> , <u><a href="#">NP_006312</a></u>   |
| Synonyms:                             | ARI2; TRIAD1  |
| Storage:                              | Room Temperature, or 2-8°C for long term storage  |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <u><a href="#">NP_006312</a></u>  |
| Locus ID:                             | 10425   |
| Cytogenetics:                         | 3p21.31   |
| Protein Families:                     | Druggable Genome  |



[View online »](#)