

Product datasheet for **TB401933**

PIM2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PIM2 (NM_006875) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201933
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	PIM2 Mouse Monoclonal Antibody [Clone ID: OTI9A5] (TA501166)
ACCN:	NM_006875 , NP_006866
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006866
Locus ID:	11040
Cytogenetics:	Xp11.23
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Acute myeloid leukemia



[View online »](#)