

Product datasheet for **TB402013**

MTGR1 (CBFA2T2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CBFA2T2 (NM_001039709), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202013
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001039709 , NP_001034798
Synonyms:	EHT; MTGR1; p85; ZMYND3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001034798
Locus ID:	9139
Cytogenetics:	20q11.21-q11.22
Protein Families:	Transcription Factors



[View online »](#)