

Product datasheet for **TB402076**

MLC1SA (MYL6B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MYL6B (NM_002475), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202076
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MLC1SA (MYL6B) Mouse Monoclonal Antibody [Clone ID: OTI7F10] (TA811146)
ACCN:	NM_002475 , NP_002466
Synonyms:	MLC1SA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002466
Locus ID:	140465
Cytogenetics:	12q13.2
Protein Pathways:	Vascular smooth muscle contraction



[View online »](#)