

Product datasheet for **TB402116**

TM7SF2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TM7SF2 (NM_003273) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202116
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OT14C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_003273 , NP_003264
Synonyms:	ANG1; C14SR; DHCR14A; NET47
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003264
Locus ID:	7108
Cytogenetics:	11q13.1
Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, Steroid biosynthesis



[View online »](#)