

Product datasheet for **TB402420**

Serine Palmitoyltransferase (SPTLC2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SPTLC2 (NM_004863) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202420
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_004863 , NP_004854
Synonyms:	hLCB2a; HSN1C; LCB2; LCB2A; NSAN1C; SPT2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004854
Locus ID:	9517
Cytogenetics:	14q24.3
Protein Pathways:	Metabolic pathways, Sphingolipid metabolism



[View online »](#)