

Product datasheet for **TB402438**

KLF5 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KLF5 (NM_001730) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202438
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	KLF5 Mouse Monoclonal Antibody [Clone ID: OTI6D2] (TA811868)
ACCN:	<u>NM_001730, NP_001721</u>
Synonyms:	BTEB2; CKLF; IKLF
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001721</u>
Locus ID:	688
Cytogenetics:	13q22.1
Protein Families:	Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transcription Factors



[View online »](#)