

Product datasheet for **TB402474**

TLE2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TLE2 (NM_003260), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202474
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TLE2 Mouse Monoclonal Antibody [Clone ID: OTI5A6] (TA504240)
ACCN:	<u>NM_003260</u> , <u>NP_003251</u>
Synonyms:	ESG; ESG2; GRG2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_003251</u>
Locus ID:	7089
Cytogenetics:	19p13.3
Protein Families:	Transcription Factors



[View online »](#)