

Product datasheet for **TB402475**

Phosphoserine Aminotransferase (PSAT1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PSAT1 (NM_058179), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202475
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_058179 , NP_478059
Synonyms:	EPIP; NLS2; PSA; PSAT; PSATD
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_478059
Locus ID:	29968
Cytogenetics:	9q21.2
Protein Pathways:	Glycine, serine and threonine metabolism, Metabolic pathways, Vitamin B6 metabolism



[View online »](#)