

## Product datasheet for **TB402506**

### MCM2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MCM2 (NM_004526) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202506
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	MCM2 Mouse Monoclonal Antibody [Clone ID: OTI8A11] (TA805566)
ACCN:	<a href="#">NM_004526</a> , <a href="#">NP_004517</a>
Synonyms:	BM28; CCNL1; cdc19; CDCL1; D3S3194; DFNA70; MITOTIN
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_004517</a>
Locus ID:	4171
Cytogenetics:	3q21.3
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transcription Factors
Protein Pathways:	Cell cycle, DNA replication



[View online »](#)