

Product datasheet for **TB402720**

Carboxypeptidase A (CPA1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CPA1 (NM_001868) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202720
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Carboxypeptidase A (CPA1) Mouse Monoclonal Antibody [Clone ID: OTI2A3] (TA500053)
ACCN:	NM_001868 , NP_001859
Synonyms:	CPA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001859
Locus ID:	1357
Cytogenetics:	7q32.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Protease, Secreted Protein



[View online »](#)