

Product datasheet for **TB402847**

DDIT4 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DDIT4 (NM_019058) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202847
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_019058 , NP_061931
Synonyms:	Dig2; REDD-1; REDD1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_061931
Locus ID:	54541
Cytogenetics:	10q22.1
Protein Pathways:	mTOR signaling pathway



[View online »](#)